

Customer Notice

Moog Medical Devices Group
4314 Zevex Park Lane
Salt Lake City, Utah 84123

ph 801.264.1001
fx 801.264.1051
800.970.2337

www.moog.com/medical

Thank you for your purchase of the Moog Infinity Set. For customers that have purchased an enteral set from Moog before, we have made some modifications to ensure they meet and surpass your needs as an Infinity user. Below is a diagram outlining improved features:

Backpack Strap Loops


The loops offer ease in securing the set during ambulatory use. When using the set within a backpack, thread the Velcro® straps through the loops, then around the funnel and secure tightly.

Cap & Funnel

The easy-to-secure cap and funnel is the most noticeable change to the set. Unlike other caps, the new design requires a quarter turn to ensure a secure closure.

Anti-Roll Flange

The addition of this new feature reduces the chance of the bag folding over on the tubing exit port.



It is Moog's intention to transition to the new sets as quickly as possible. For a period of time, a possibility does exist that you will receive a case of sets with the former design. The performance of both sets, regardless of which one you receive, is consistent with Moog's commitment to our customers for quality, flexible and easy to use products.

If you have any questions, please contact our Technical Support Team at (800) 970-2337. Thank you for being a Moog customer and we look forward to serving you in the future.