

Life Science Laboratories - Rush

7500 West Henrietta Road | Rush, NY 14543 | 585.533.1672 | fax 585.533.1796

Send Final Report To:

Company _____
 Attention _____
 Address _____
 City _____ State _____
 Zip _____ Country _____
 Phone _____ Fax _____
 Email _____

Billing Information:

Address (If different than report to) _____

 Purchase Order # _____
 Quote # / Cost _____
 Sample Replacement Value _____

Sample Information

Test Name/Guideline (ISO, AAMI, USP, etc.)	Number of Samples	Product Name	Lot Number	Other Identifier

Storage Conditions

- Room temperature
- Refrigerate (2-8°C)
- Freeze (-10 to -25°C)
- Other

Submitted Material Type

- Plastic/polymer/metal
- Elastomer
- Coated material, composite, laminate
- Pharmaceuticals
- Other

Safety Precautions

- None/unknown (standard precautions will be used)
- MSDS enclosed
- Flammable
- Biohazard level
- Other

Sample Disposition*

- Discard
- Return unused sample
- Return used and unused sample

*Product will be destroyed unless otherwise indicated

Return Instructions (if applicable)

- UPS
 - FedEx
 - Other
- Account # _____

Test Requirements (please complete where applicable)

Extraction Condition Options

- 121°C - 1 hour
- 70°C - 24 hour
- 50°C - 72 hour
- 37°C - 72 hour (not applicable for USP Class plastics)
- 37°C - 24 hour (cytotoxicity testing in SS MEM)
- Other:

Extracting Media Options

- Saline & vegetable oil (ISO)
- Saline, vegetable oil, 1:20 alcohol:saline, & polyethylene glycol (Class VI USP plastics only)
- Serum supplemented MEM (Elution-cytotoxicity testing)
- Other:

Sensitization/Irritation Evaluation

- Direct contact
- Extract

Do NOT sub-divide (cut) test sample for extraction (Test sample will be cut into the appropriate sizes needed for testing unless box is checked.)

Do NOT test entire article (Entire sample will be tested unless box is checked. Identify specific components or materials to be excluded below.):

Samples Submitted are:

- Sterile, processed by indicated method → → →
- Non-sterile
 - process per selected method → → → →
 - not intended to be sold or used sterile for clinical use*

Sterilization Method

- EtO
- steam
- gamma
- other

Implantation Evaluation*

- histopathology (ISO)
- macroscopic

NOTE*: Samples for histology and post mortem clinical evaluation will be submitted to the Division of Laboratory Animal Medicine - University of Rochester. Samples that are sold and used sterile shall be submitted sterile or be sterilized by the intended method prior to testing. Samples submitted for implantation testing must be submitted sterile or be sterilized before testing.

Results to be Released By

- email (pdf)
- fax
- US mail
- UPS or FedEx account:

Sponsor contact for questions pertaining to testing (print name): _____

Comments or Special Instructions _____

Sponsor Authorization _____ Date _____

Unless other arrangements have been made, payment terms are net 30, quotes are valid for 30 days, prices are FOB Moog Medical Devices Group, Life Sciences Laboratories and are valid for 6 months after acceptance of the quote.

Internal Use Only	
Received By _____	Date _____
Approved By _____	Date _____
Test Number(s): _____	