


Pressure Infusion

PRODUCT NAME	DESCRIPTION	CATALOG	UOM
Infu-Surg® Disposable Pressure Infusion Bags	500ml Pressure Infuser Bag	4005	5 / Box 25 / Case
	500ml Pressure Infuser Bag With Hook	4005H	5 / Box 25 / Case
	1000ml Pressure Infuser Bag	4010	5 / Box 25 / Case
	1000ml Pressure Infuser Bag With Hook	4010H	5 / Box 25 / Case
	3000ml Pressure Infuser Bag	4030	5 / Box
	Surgi-Press®	500cc Reusable Pressure Infuser	4081



Surgical Supplies

PRODUCT NAME	DESCRIPTION	CATALOG	UOM	
Surgi-Kit® Instrument Holders and Suture Organizers	8 Pocket without tape: Sterile - 19cm x 23cm	SK-100	50 / Box 200 / Case	
	8 Pocket without tape: Non-Sterile - 19cm x 23cm	SK-100N	500 / Case	
	4 Pocket without tape: Sterile - 19cm x 11cm	SK-150	100 / Box 400 / Case	
	4 Pocket without tape: Non-Sterile - 19cm x 11cm	SK-150N	500 / Case	
	2 Pocket without tape: Sterile - 28cm x 18cm	SK-200	50 / Box 200 / Case	
	2 Pocket without tape: Non-Sterile - 28cm x 18cm	SK-200N	500 / Case	
	8 Pocket with tape: Sterile - 19cm x 23cm	SK-300	50 / Box 200 / Case	
	8 Pocket with tape: Non-Sterile - 19cm x 23cm	SK-300N	500 / Case	
	4 Pocket with tape: Sterile - 19cm x 11cm	SK-350	75 / Box 300 / Case	
	4 Pocket with tape: Non-Sterile - 19cm x 11cm	SK-350N	500 / Case	
	2 Pocket with tape: Sterile - 28cm x 18cm	SK-400	50 / Box 200 / Case	
	2 Pocket with tape: Non-Sterile - 28cm x 18cm	SK-400N	500 / Case	
	10 Pocket with tape: Sterile - 15cm x 19cm	SK-650	50 / Box 200 / Case	
	10 Pocket with tape: Non-Sterile - 15cm x 19cm	SK-650N	500 / Case	
	Lap-Kit® Laparoscopic Instrument Holders	2 Pocket with tape: Sterile - 68mm x 30cm	LK-700	25 / Box 100 / Case
		2 Pocket with tape: Non-Sterile - 68mm x 30cm	LK-700N	200 / Case
4 Pocket with tape: Sterile - 33mm x 30cm		LK-800	25 / Box 100 / Case	
3 Pocket with tape: Sterile - (1)68mm x 30cm (2)33mm x 30cm		LK-900	25 / Box 100 / Case	
3 Pocket with tape: Non-Sterile - (1)68mm x 30cm (2)33mm x 30cm		LK-900N	200 / Case	



Reusable Blood Pressure Cuffs

PRODUCT NAME

**Surgi-Cuf®
Reusable Blood
Pressure Cuffs**



DESCRIPTION

CATALOG

UOM

Neonatal Size (8-14cm)

1-Tube Female Luer Lock Connector	7105	5 / Box
1-Tube Multi-purpose Screw Connector	7305	5 / Box
1-Tube Rectus Style Quick Disconnect	7405	5 / Box
1-Tube Female CPC Connector	7605	5 / Box
1-Tube Male CPC Connector	7705	5 / Box
2-Tube Multipurpose Screw Connector	7205	5 / Box
2-Tube Male and Female CPC Connector	7805	5 / Box

Child Size (13-20cm)

1-Tube Female Luer Lock Connector	7106	5 / Box
1-Tube Multi-purpose Screw Connector	7306	5 / Box
1-Tube Rectus Style Quick Disconnect	7406	5 / Box
1-Tube Female CPC Connector	7606	5 / Box
1-Tube Male CPC Connector	7706	5 / Box
2-Tube Multipurpose Screw Connector	7206	5 / Box
2-Tube Male and Female CPC Connector	7806	5 / Box

Small Adult Size (18-26cm)

1-Tube Female Luer Lock Connector	7107	5 / Box
1-Tube Multi-purpose Screw Connector	7307	5 / Box
1-Tube Rectus Style Quick Disconnect	7407	5 / Box
1-Tube Female CPC Connector	7607	5 / Box
1-Tube Male CPC Connector	7707	5 / Box
2-Tube Multipurpose Screw Connector	7207	5 / Box
2-Tube Male and Female CPC Connector	7807	5 / Box

Adult Size (26-35cm)

1-Tube Female Luer Lock Connector	7108	5 / Box
1-Tube Multi-purpose Screw Connector	7308	5 / Box
1-Tube Rectus Style Quick Disconnect	7408	5 / Box
1-Tube Female CPC Connector	7608	5 / Box
1-Tube Male CPC Connector	7708	5 / Box
2-Tube Multipurpose Screw Connector	7208	5 / Box
2-Tube Male and Female CPC Connector	7808	5 / Box

Adult Plus Size (32-42cm)

1-Tube Female Luer Lock Connector	7108P	5 / Box
1-Tube Multi-purpose Screw Connector	7308P	5 / Box
1-Tube Rectus Style Quick Disconnect	7408P	5 / Box
1-Tube Female CPC Connector	7608P	5 / Box
1-Tube Male CPC Connector	7708P	5 / Box
2-Tube Multipurpose Screw Connector	7208P	5 / Box
2-Tube Male and Female CPC Connector	7808P	5 / Box





Large Adult Size (32-42cm)

1-Tube Female Luer Lock Connector	7109	5 / Box
1-Tube Multi-purpose Screw Connector	7309	5 / Box
1-Tube Rectus Style Quick Disconnect	7409	5 / Box
1-Tube Female CPC Connector	7609	5 / Box
1-Tube Male CPC Connector	7709	5 / Box
2-Tube Multipurpose Screw Connector	7209	5 / Box
2-Tube Male and Female CPC Connector	7809	5 / Box

Thigh Model (42-50cm)

1-Tube Female Luer Slip Connector	7110	5 / Box
1-Tube Multi-purpose Screw Connector	7310	5 / Box
1-Tube Rectus Style Quick Disconnect	7410	5 / Box
1-Tube Female CPC Connector	7610	5 / Box
1-Tube Male CPC Connector	7710	5 / Box
2-Tube Multipurpose Screw Connector	7210	5 / Box
2-Tube Male and Female CPC Connector	7810	5 / Box

Vinyl, Neo-Natal and Soft Disposable Blood Pressure Cuffs

PRODUCT NAME	DESCRIPTION	CATALOG	UOM	
Surgi-Cuf® Adult and Child Vinyl Blood Pressure Cuffs 	Child Size (15-22cm)			
	1-Tube Female Luer Lock Connector	5106	10 / Box	
	1-Tube Quick-disconnect Connector	5406	10 / Box	
	2-Tube Multi-purpose Screw Connector	5206	10 / Box	
	Small Adult Size (17-25cm)			
	1-Tube Female Luer Lock Connector	5107	10 / Box	
	1-Tube Quick-disconnect Connector	5407	10 / Box	
	2-Tube Multi-purpose Screw Connector	5207	10 / Box	
	Adult Size (24-32cm)			
	1-Tube Female Luer Lock Connector	5108	10 / Box	
	1-Tube Quick-disconnect Connector	5408	10 / Box	
	2-Tube Multi-purpose Screw Connector	5208	10 / Box	
	Adult Plus Size (28-37cm)			
	1-Tube Female Luer Lock Connector	5108P	10 / Box	
	1-Tube Quick-disconnect Connector	5408P	10 / Box	
	2-Tube Multi-purpose Screw Connector	5208P	10 / Box	
	Large Adult Size (32-42cm)			
	1-Tube Female Luer Lock Connector	5109	10 / Box	
	1-Tube Quick-disconnect Connector	5409	10 / Box	
	2-Tube Multi-purpose Screw Connector	5209	10 / Box	
Thigh Model (42-50cm)				
1-Tube Female Luer Lock Connector	5110	10 / Box		
1-Tube Quick-disconnect Connector	5410	10 / Box		
2-Tube Multi-purpose Screw Connector	5210	10 / Box		
Surgi-Cuf® Sterile Neonatal Disposable Blood Pressure Cuffs 	Size 1 (3-6cm)			
	1-Tube Male Luer Slip Connector	5101	10 / Box	
	2-Tube Male Luer Slip Connector	5201	10 / Box	
	Size 2 (6-9.5cm)			
	1-Tube Male Luer Slip Connector	5102	10 / Box	
	2-Tube Male Luer Slip Connector	5202	10 / Box	
	Size 3 (8-12cm)			
	1-Tube Male Luer Slip Connector	5103	10 / Box	
	2-Tube Male Luer Slip Connector	5203	10 / Box	
	Size 4 (9-14cm)			
	1-Tube Male Luer Slip Connector	5104	10 / Box	
	2-Tube Male Luer Slip Connector	5204	10 / Box	
	Size 5 (12-18cm)			
	1-Tube Male Luer Slip Connector	5105	10 / Box	
	2-Tube Male Luer Slip Connector	5205	10 / Box	
Surgi-Cuf® Soft Disposable Blood Pressure Cuffs 	Child (13-20cm)	1-Tube Quick-disconnect Connector	6406	10 / Box
	Small Adult (18-26cm)	1-Tube Quick-disconnect Connector	6407	10 / Box
	Adult (26-35cm)	1-Tube Quick-disconnect Connector	6408	10 / Box
	Adult Long (29-38cm)	1-Tube Quick-disconnect Connector	6408P	10 / Box
	Large Adult (32-42cm)	1-Tube Quick-disconnect Connector	6409	10 / Box
	Thigh (42-50cm)	1-Tube Quick-disconnect Connector	6410	10 / Box
	Surgi-Cuf® Blood Pressure Adapter Standardization System 	Inflation Bulbs w/ Male Luer	5000B	5 / Box
1-Tube Adapters w/ Male Quick Connect + Male Luer Lock (HP Monitors)		5002	10 / Box	
2-Tube Adapters w/ Male Luer and Female CPC (Marquette Monitors)		5004	10 / Box	
2-Tube Adapters w/ Male luer and Male + Female CPC (Datex Monitors)		5006	10 / Box	
Female Quick Connect to Male and Female CPC		5007	5 / Pack	
Female Quick Connect to Two Multi Purpose Screws		5008	5 / Pack	
Female Quick Connect to Two Female CPC		5009	5 / Pack	
Female Quick Connect to Two 1/8" Tubing w/Bulb + Bare Tube and barb		5010	5 / Pack	
Female Quick Connect to Two 1/8" Tubing w/ Barbs		5011	5 / Pack	
Male Quick Connect to Two Tube w/ Female Luer Lock		5012	5 / Pack	
Female Quick Connect to Two 1/8" Tubes w/ Bulb + Bare Tube		5018	5 / Pack	
Female Quick Connect Adapter for 1/8" Tubing		5050	5 / Pack	
Female Quick Connect Adapter for 3/16" Tubing		5051	5 / Pack	
Female Quick Connect to Female Luer Lock		5052	5 / Pack	
Female Quick Connect to Multi-Purpose Screw		5053	5 / Pack	
Female Quick Connect to Female CPC		5054	5 / Pack	
Female Quick Connect to Male CPC		5055	5 / Pack	
Female Quick Connect to Male Nipset		5056	5 / Pack	
Male Quick Connect to Female Luer Lock		5057	5 / Pack	

Gastric Lavage

PRODUCT NAME	DESCRIPTION	CATALOG	UOM
Tum-E-Vac® Gastric Lavage Kits with Lavage Tubes and/or Charcoal	Lavage Kit with 3L bag	2055	10 / Box
	Lavage Kit with spike, no bag with tubes and / or charcoal	2075	10 / Box
	32fr Lavage Tubes - Lavage Kit w/32fr & 3L bag	2054	10 / Box
	Lavage Kit w/32fr Lavage Tubes & Aqua & 2-25g Activated Charcoal	2057	10 / Box
	Lavage Kit w/32fr Lavage Tubes & (2) 25g Activated Charcoal and (2) 25g Activated Charcoal w/Sorbitol	2058	10 / Box
	Lavage Kit w/32fr Lavage Tubes & (1) 25g Activated Charcoal and (1) 25g Activated Charcoal w/Sorbitol	2059	10 / Box
	Lavage Kit w/ Spike, Sorbitol & Aqua and (2) 25g Activated Charcoal w/Sorbitol	2081	10 / Box



Tum-E-Vac® Gastric Lavage Tubes	24fr Lavage Tubes	2124	10 / Box
	28fr Lavage Tubes	2128	10 / Box
	32fr Lavage Tubes	2132	10 / Box
	36fr Lavage Tubes	2136	10 / Box



Charcoal and Bite Block	Bite Block	2060	20 / Box
	Actidose-AQUA	3010	12 / Box
	Actidose-SORBITOL	3020	12 / Box



Wound Irrigation

Combiport® & Combiguard II® Wound Irrigation

Combiport®/Combiguard® Kit	10100	50 / Box
Combiguard® II - use with 60cc syringe	10222	50 / Box
Comiguard II® - Non-Sterile	10222NS	500 / Box



Moog Medical Devices Group

251 Seneca Street Buffalo, NY 14204
1.800.521.1022 716.842.4000 fax 716.842.2546
med_hp@moog.com www.moog.com

MOOG